

Product datasheet for **TA376113**

TENT5C Rabbit Polyclonal Antibody


Product data:

| | |
|-------------------------|--|
| Product Type: | Primary Antibodies |
| Applications: | WB |
| Recommended Dilution: | WB,1:500 - 1:2000 |
| Reactivity: | Human |
| Modifications: | Unmodified |
| Host: | Rabbit |
| Isotype: | IgG |
| Clonality: | Polyclonal |
| Immunogen: | Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human FAM46C (NP_060179.2). |
| Formulation: | Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| Concentration: | lot specific |
| Purification: | Affinity purification |
| Conjugation: | Unconjugated |
| Storage: | Store at -20°C. Avoid freeze / thaw cycles. |
| Stability: | Shelf life: one year from despatch. |
| Predicted Protein Size: | 45kDa |
| Gene Name: | family with sequence similarity 46 member C |
| Database Link: | Entrez Gene 54855 Human Q5VWP2 |
| Background: | Nucleotidyltransferase that act as a non-canonical poly(A RNA polymerase which enhances mRNA stability and gene expression. Mainly targets mRNAs encoding endoplasmic reticulum-targeted protein and may be involved in induction of cell death. |
| Synonyms: | FLJ20202 |



[View online »](#)

Product images:



Western blot analysis of extracts of various cell lines, using FAM46^C antibody (TA376113) at 1:1000 dilution. | Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. | Lysates/proteins: 25ug per lane. | Blocking buffer: 3% nonfat dry milk in TBST. | Detection: ECL Basic Kit . | Exposure time: 30s.