

## Product datasheet for **TA376117S**


### CIAO2B Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	WB, 1:500 - 1:2000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Reactivity:	Human, Mouse
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	18kDa
Gene Name:	family with sequence similarity 96 member B
Database Link:	<a href="#">Entrez Gene 51647 Human</a> <a href="#">Q9Y3D0</a>
Background:	Involved in chromosome segregation; iron-sulfur cluster assembly; and protein maturation by iron-sulfur cluster transfer. Located in cytosol; nucleoplasm; and spindle. Part of CIA complex and MMXD complex.
Synonyms:	CGI-128



[View online »](#)

**Product images:**


Western blot analysis of various lysates using FAM96B Rabbit pAb (TA376117) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 15s.