

## Product datasheet for **TA376136S**

### **FASTKD5 Rabbit Polyclonal Antibody**


#### **Product data:**

Product Type:	Primary Antibodies
Applications:	ELISA, ICC/IF, WB
Recommended Dilution:	WB, 1:500 - 1:2000 IF/ICC, 1:50 - 1:200 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Reactivity:	Human, Mouse, Rat
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	87kDa
Gene Name:	FAST kinase domains 5
Database Link:	<a href="#">Entrez Gene 60493 Human Q7L8L6</a>
Background:	Enables rRNA binding activity. Involved in mitochondrial RNA processing. Located in mitochondrial nucleoid and ribonucleoprotein granule.
Synonyms:	dj1187M17.5; FLJ13149; FLJ58294; KIAA1792



[View online »](#)

## Product images:



Western blot analysis of various lysates using FASTKD5 Rabbit pAb (TA376136) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 90s.