

Product datasheet for **TA376308**

FNDC3A Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:200 - 1:2000
Reactivity:	Human, Mouse
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 900-1170 of human FNDC3A (NP_001265367.1).
Formulation:	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	125kDa/131kDa
Gene Name:	fibronectin type III domain containing 3A
Database Link:	Entrez Gene 22862 Human Q9Y2H6
Background:	Mediates spermatid-Sertoli adhesion during spermatogenesis.
Synonyms:	bA203I16.1; bA203I16.5; FLJ31509; FNDC3; HUGO; KIAA0970



[View online »](#)