

## Product datasheet for **TA376370S**

### FTO Rabbit Polyclonal Antibody

#### Product data:


Product Type:	Primary Antibodies
Applications:	IP, WB
Recommended Dilution:	WB,1:500 - 1:2000 IP,1:50 - 1:200
Reactivity:	Human, Mouse, Rat
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 416-505 of human FTO (NP_001073901.1).
Formulation:	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	6kDa/12kDa/14kDa/58kDa
Gene Name:	fat mass and obesity associated
Database Link:	<a href="#">Entrez Gene 79068 Human Q9C0B1</a>
Background:	This gene is a nuclear protein of the AlkB related non-haem iron and 2-oxoglutarate-dependent oxygenase superfamily but the exact physiological function of this gene is not known. Other non-heme iron enzymes function to reverse alkylated DNA and RNA damage by oxidative demethylation. Studies in mice and humans indicate a role in nervous and cardiovascular systems and a strong association with body mass index, obesity risk, and type 2 diabetes.




[View online »](#)

Synonyms: KIAA1752; MGC5149

**Product images:**



Western blot analysis of extracts of various cell lines, using FTO antibody ([TA376370]) at 1:1000 dilution. | Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. | Lysates/proteins: 25ug per lane. | Blocking buffer: 3% nonfat dry milk in TBST. | Detection: ECL Basic Kit . | Exposure time: 30s.



Immunoprecipitation analysis of 300ug extracts of SH-SY5Y cells using 3ug FTO antibody ([TA376370]). Western blot was performed from the immunoprecipitate using FTO antibody ([TA376370]) at a dilution of 1:1000.