

Product datasheet for **TA376667**


GNPNAT1 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	ICC/IF, WB
Recommended Dilution:	WB,1:500 - 1:2000 IF,1:50 - 1:200
Reactivity:	Human, Mouse, Rat
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-91 of human GNPAT1 (NP_932332.1).
Formulation:	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Gene Name:	glucosamine-phosphate N-acetyltransferase 1
Database Link:	Entrez Gene 64841 Human Q96EK6
Synonyms:	FLJ10607; GNA1; GNPAT; Gpnat1



[View online »](#)

Product images:


Western blot analysis of extracts of various cell lines, using GNPAT1 antibody (TA376667) at 1:1000 dilution. | Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. | Lysates/proteins: 25ug per lane. | Blocking buffer: 3% nonfat dry milk in TBST. | Detection: ECL Basic Kit . | Exposure time: 15s.



Immunofluorescence analysis of L929 cells using GNPAT1 antibody (TA376667) at dilution of 1:100. Blue: DAPI for nuclear staining.