

Product datasheet for **TA376901**

H3FA (HIST1H3A) Rabbit Polyclonal Antibody

Product data:

| | |
|-------------------------|--|
| Product Type: | Primary Antibodies |
| Applications: | Dot, WB |
| Recommended Dilution: | DB, 1:500 - 1:4000 WB, 1:500 - 1:2000 |
| Reactivity: | Human, Mouse, Rat, Other (Wide Range) |
| Host: | Rabbit |
| Isotype: | IgG |
| Clonality: | Polyclonal |
| Immunogen: | Recombinant protein of human Lactic acid-Histone H3-K18. |
| Formulation: | Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3. |
| Concentration: | lot specific |
| Purification: | Affinity purification |
| Conjugation: | Unconjugated |
| Storage: | Store at -20°C. Avoid freeze / thaw cycles. |
| Stability: | Shelf life: one year from despatch. |
| Predicted Protein Size: | 15kDa |
| Gene Name: | histone cluster 1, H3a |
| Database Link: | Entrez Gene 8350 Human P68431 |


Background: Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. This structure consists of approximately 146 bp of DNA wrapped around a nucleosome, an octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H3 family. Transcripts from this gene lack polyA tails; instead, they contain a palindromic termination element. This gene is found in the large histone gene cluster on chromosome 6p22-p21.3.




[View online »](#)

Synonyms: H3/A; H3/b; H3/c; H3/d; H3/f; H3/h; H3/i; H3/j; H3/k; H3/l; H3FA; H3FB; H3FC; H3FD; H3FF; H3FH; H3FI; H3FJ; H3FK; H3FL

Product images:



Dot-blot analysis of all sorts of peptides using Lactic acid-Histone H3-K18 antibody (TA376901) at 1:1000 dilution.



Western blot analysis of extracts of MCF7 cells, using Lactic acid-Histone H3-K18 antibody (TA376901) at 1:500 dilution. MCF7 cells were treated by Sodium lactate (100mM) for 24h. | Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. | Lysates/proteins: 25ug per lane. | Blocking buffer: 3% BSA. | Detection: ECL Basic Kit. | Exposure time: 5min.