

# **Product datasheet for TA377052**

# H2BC5 Rabbit Polyclonal Antibody

# **Product data:**

#### OriGene Technologies, Inc.


9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Primary Antibodies
Applications:	ICC/IF, IHC, IP, WB
Recommended Dilution:	WB,1:500 - 1:2000 IHC,1:50 - 1:200 IF,1:50 - 1:100 IP,1:50 - 1:200
Reactivity:	Human, Mouse, Rat, Other (Wide Range)
Modifications:	Acetyl K15
Host:	Rabbit
lsotype:	lgG
Clonality:	Polyclonal
Immunogen:	A synthetic acetylated peptide corresponding to residues surrounding K15 of human Histone H2B
Formulation:	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Gene Name:	histone cluster 1, H2bd
Database Link:	<u>Entrez Gene 3017 Human</u> <u>P62807</u>
Synonyms:	dJ221C16.6; H2B.1B; H2B.b; H2B/b; H2BFB; HIRIP2; MGC90432; OTTHUMP00000016143




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

## **Product images:**




Western blot analysis of extracts of various cell lines, using Acetyl-Histone H2B-K15 antibody (TA377052) at 1:1000 dilution.HeLa cells were treated by TSA (1 uM) at 37°C for 18 hours.NIH/3T3 cells were treated by TSA (1 uM) at 37°C for 18 hours.C6 cells were treated by TSA (1 uM) at 37°C for 18 hours.|Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.|Lysates/proteins: 25ug per lane.|Blocking buffer: 3% BSA.|Detection: ECL Basic Kit.|Exposure time: 10s.




Immunofluorescence analysis of C6 cells using Acetyl-Histone H2B-K15 antibody (TA377052) at dilution of 1:100.C6 cells were treated by TSA (1 uM) at 37°C for 18 hours. Blue: DAPI for nuclear staining.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

## 



Immunofluorescence analysis of HeLa cells using Acetyl-Histone H2B-K15 antibody (TA377052) at dilution of 1:100.HeLa cells were treated by TSA (1 uM) at 37°C for 18 hours. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of NIH/3T3 cells using Acetyl-Histone H2B-K15 antibody (TA377052) at dilution of 1:100.NIH/3T3 cells were treated by TSA (1 uM) at 37°C for 18 hours. Blue: DAPI for nuclear staining.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US