

Product datasheet for **TA377301**


IAH1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	WB, 1:500 - 1:1000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Reactivity:	Human, Mouse, Rat
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	28kDa
Gene Name:	isoamyl acetate-hydrolyzing esterase 1 homolog
Database Link:	Entrez Gene 285148 Human Q2TAA2
Background:	Enables identical protein binding activity. Predicted to be involved in lipid catabolic process.
Synonyms:	MGC102860



[View online »](#)

Product images:


Western blot analysis of various lysates