

Product datasheet for **TA377334S**


IFI35 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	WB, 1:500 - 1:2000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Reactivity:	Human, Mouse
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	32kDa
Gene Name:	interferon induced protein 35
Database Link:	Entrez Gene 3430 Human P80217
Background:	Enables identical protein binding activity. Involved in several processes, including macrophage activation involved in immune response; positive regulation of defense response; and regulation of signal transduction. Located in several cellular components, including cytosol; extracellular space; and nucleus.
Synonyms:	FLJ21753; IFP35



[View online »](#)

Product images:

Western blot analysis of various lysates using IFI35 Rabbit pAb ([TA377334]) at 1:1000 dilution. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.