

Product datasheet for **TA378047**


LIMD1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	WB, 1:500 - 1:2000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Reactivity:	Human, Mouse
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	72kDa
Gene Name:	LIM domains containing 1
Database Link:	Entrez Gene 8994 Human Q9UGP4
Background:	Predicted to enable transcription corepressor activity. Involved in several processes, including negative regulation of hippo signaling; regulation of gene expression; and response to hypoxia. Acts upstream of or within P-body assembly and gene silencing by miRNA. Located in several cellular components, including P-body; adherens junction; and focal adhesion. Part of RISC complex.
Synonyms:	LIMD1



[View online »](#)

Product images:


Western blot analysis of various lysates using LIMD1 Rabbit pAb (TA378047) at 1:1000 dilution. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 10s.