

Product datasheet for **TA378061**


Lins1 Rabbit Polyclonal Antibody

Product data:

| | |
|-------------------------|--|
| Product Type: | Primary Antibodies |
| Applications: | ELISA, WB |
| Recommended Dilution: | WB, 1:500 - 1:2000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. |
| Reactivity: | Mouse, Rat |
| Modifications: | Unmodified |
| Host: | Rabbit |
| Isotype: | IgG |
| Clonality: | Polyclonal |
| Formulation: | Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3. |
| Concentration: | lot specific |
| Purification: | Affinity purification |
| Conjugation: | Unconjugated |
| Storage: | Store at -20°C. Avoid freeze / thaw cycles. |
| Stability: | Shelf life: one year from despatch. |
| Predicted Protein Size: | 55kDa/85kDa |
| Gene Name: | lines homolog 1 |
| Database Link: | Entrez Gene 72635 Mouse Q3U1D0 |
| Background: | Predicted to be involved in cognition. Human ortholog(s) of this gene implicated in autosomal recessive non-syndromic intellectual disability. Orthologous to human LINS1 (lines homolog 1). |
| Synonyms: | FLJ10583; WINS1 |



[View online »](#)

Product images:


Western blot analysis of various lysates using LINS1 Rabbit pAb (TA378061) at 1:1000 dilution. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 1s.