

Product datasheet for **TA378567S**


AF9 (MLLT3) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	WB, 1:500 - 1:2000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Reactivity:	Human
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	63kDa
Gene Name:	myeloid/lymphoid or mixed-lineage leukemia; translocated to, 3
Database Link:	Entrez Gene 4300 Human P42568
Background:	Enables chromatin binding activity and lysine-acetylated histone binding activity. Involved in several processes, including hematopoietic stem cell differentiation; positive regulation of transcription, DNA-templated; and regulation of stem cell division. Acts upstream of or within negative regulation of canonical Wnt signaling pathway and positive regulation of Wnt signaling pathway, planar cell polarity pathway. Located in cytosol and nucleoplasm. Part of transcription elongation factor complex.
Synonyms:	AF-9; AF9; FLJ2035; YEATS3



[View online »](#)

Product images:


Western blot analysis of various lysates using MLLT3/AF9 Rabbit pAb ([TA378567]) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 90s.