

Product datasheet for **TA378623**

MPRIP Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:200 - 1:2000
Reactivity:	Human
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-310 of human MPRIP (NP_055949.2).
Formulation:	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	114kDa/116kDa/118kDa
Gene Name:	myosin phosphatase Rho interacting protein
Database Link:	Entrez Gene 23164 Human Q6WCQ1
Background:	MPRIP, also named as KIAA0864, MRIP, RHOIP3 and p116Rip, is required for the regulation of the actin cytoskeleton by RhoA and ROCK1. MPRIP targets myosin phosphatase to the actin cytoskeleton. Depletion leads to an increased number of stress fibers in smooth muscle cells through stabilization of actin fibers by phosphorylated myosin. Overexpression of MRIP as well as its F-actin-binding region leads to disassembly of stress fibers in neuronal cells. The antibody recognizes the C-term of isoform2 and isoform3 of MPRIP.
Synonyms:	KIAA0864; M-RIP; MRIP; p116Rip; RHOIP3; RIP3



[View online »](#)