


## Product datasheet for TA379194S

### NOC4L Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	WB, 1:200 - 1:2000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Reactivity:	Human
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	58kDa
Gene Name:	nucleolar complex associated 4 homolog
Database Link:	<a href="#">Entrez Gene 79050 Human</a> <a href="#">Q9BVI4</a>
Background:	Enables RNA binding activity. Predicted to be involved in rRNA processing. Located in nucleolus and nucleoplasm.
Synonyms:	MGC3162; NET49; NOC4


[View online »](#)

**Product images:**


Western blot analysis of various lysates using NOC4L Rabbit pAb (TA379194) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 90s.