

Product datasheet for **TA379317S**


NRN1L Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:500 - 1:2000
Reactivity:	Human, Mouse, Rat
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant protein of human NRN1L
Formulation:	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	17kDa
Gene Name:	neuritin 1 like
Database Link:	Entrez Gene 123904 Human Q496H8
Background:	The protein encoded by this gene is extracellular and enhances both neurite growth and neuronal survival. The encoded protein is found both as a GPI anchored membrane-bound form and as a secreted form. This activity-related ligand functions as a homodimer or heterodimer.
Synonyms:	MGC118990; MGC118993; MRCC2446; UNQ2446



[View online »](#)

Product images:

Western blot analysis of extracts of various cell lines, using NRN1L antibody ([TA379317]) at 1:1000 dilution. | Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. | Lysates/proteins: 25ug per lane. | Blocking buffer: 3% nonfat dry milk in TBST. | Detection: ECL Basic Kit . | Exposure time: 10s.