

Product datasheet for **TA379624**

PAQR5 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:500 - 1:2000
Reactivity:	Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant protein of human PAQR5.
Formulation:	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Gene Name:	progesterin and adipoQ receptor family member 5
Database Link:	Q9NXX6
Background:	Plasma membrane progesterone (P4 receptor coupled to G proteins. Seems to act through a G(i mediated pathway. May be involved in oocyte maturation.
Synonyms:	FLJ20190; MPRG; OTTHUMP00000164447



[View online »](#)