

## Product datasheet for **TA379988**

### PKNOX1 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:200 - 1:2000
Reactivity:	Human, Mouse
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 322-436 of human PKNOX1 (NP_004562.2).
Formulation:	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	44kDa/47kDa
Gene Name:	PBX/knotted 1 homeobox 1
Database Link:	<a href="#">Entrez Gene 5316 Human P55347</a>
Background:	Activates transcription in the presence of PBX1A and HOXA1.
Synonyms:	pknox1c; PREP1



[View online »](#)

## Product images:



Western blot analysis of extracts of 293T cells, using PKNOX1 antibody (TA379988) at 1:1000 dilution. | Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. | Lysates/proteins: 25ug per lane. | Blocking buffer: 3% nonfat dry milk in TBST. | Detection: ECL Basic Kit . | Exposure time: 30s.