

## Product datasheet for **TA380990**

### RNF185 Rabbit Polyclonal Antibody

#### Product data:


Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:500 - 1:2000
Reactivity:	Human, Mouse, Rat
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human RNF185 (NP_689480.2).
Formulation:	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	14kDa/20kDa
Gene Name:	ring finger protein 185
Database Link:	<a href="#">Entrez Gene 91445 Human Q96GF1</a>
Background:	E3 ubiquitin-protein ligase that regulates selective mitochondrial autophagy by mediating 'Lys-63'-linked polyubiquitination of BNIP1. Acts in the endoplasmic reticulum (ER-associated degradation (ERAD pathway, which targets misfolded proteins that accumulate in the endoplasmic reticulum (ER for ubiquitination and subsequent proteasome-mediated degradation. Protects cells from ER stress-induced apoptosis. Responsible for the cotranslational ubiquitination and degradation of CFTR in the ERAD pathway. Preferentially associates with the E2 enzymes UBE2J1 and UBE2J2.



[View online »](#)

**Synonyms:** BSK65-MONO1; BSK65-MONO2; BSK65-PANC1; BSK65-PANC2; BSK65-TEST1; BSK65-TEST2; BSK65-TEST3; FLJ38628

### Product images:



Western blot analysis of extracts of Mouse testis, using RNF185 Rabbit pAb (TA380990) at 1:1000 dilution. | Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. | Lysates/proteins: 25ug per lane. | Blocking buffer: 3% nonfat dry milk in TBST. | Detection: ECL Basic Kit . | Exposure time: 10s.