

Product datasheet for TA381823

alpha Synuclein (SNCA) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: ICC/IF, WB

Recommended Dilution: WB,1:500 - 1:2000

IF,1:50 - 1:200

Reactivity: Human, Mouse, Rat

Modifications: Unmodified

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 61-140 of

human alpha-Synuclein (NP_000336.1).

Formulation: Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Concentration: lot specific

Purification: Affinity purification

Conjugation: Unconjugated

Storage: Store at -20°C. Avoid freeze / thaw cycles.

Stability: Shelf life: one year from despatch.

Predicted Protein Size: 11kDa/13kDa/14kDa

Gene Name: synuclein alpha

Database Link: Entrez Gene 6622 Human

P37840

Background: Alpha-synuclein is a member of the synuclein family, which also includes beta- and gamma-

synuclein. Synucleins are abundantly expressed in the brain and alpha- and beta-synuclein inhibit phospholipase D2 selectively. SNCA may serve to integrate presynaptic signaling and

membrane trafficking. Defects in SNCA have been implicated in the pathogenesis of

Parkinson disease. SNCA peptides are a major component of amyloid plaques in the brains of

patients with Alzheimer's disease. Alternatively spliced transcripts encoding different

isoforms have been identified for this gene.



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn


Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com




Synonyms:

alpha-synuclein; MGC110988; NACP; OTTHUMP00000161559; OTTHUMP00000161561; PARK1; PARK4; PD1

Product images:




Western blot analysis of extracts of various cell lines, using α -Synuclein Polyclonal Antibody (TA381823) at 1:1000 dilution. | Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. | Lysates/proteins: 25ug per lane. | Blocking buffer: 3% nonfat dry milk in TBST. | Detection: ECL Basic Kit. | Exposure time: 1s.



Western blot analysis of extracts of various cell lines, using α-Synuclein Polyclonal Antibody (TA381823) at 1:1000 dilution. | Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. | Lysates/proteins: 25ug per lane. | Blocking buffer: 3% nonfat dry milk in TBST. | Detection: ECL Basic Kit. | Exposure time: 3min.





Immunofluorescence analysis of mouse brain using α -Synuclein Rabbit pAb (TA381823) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.