

## Product datasheet for **TA381862**

### SNX32 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	ICC/IF, WB
Recommended Dilution:	WB,1:500 - 1:2000 IF,1:50 - 1:200
Reactivity:	Human, Mouse, Rat
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 144-403 of human SNX32 (NP_689973.2).
Formulation:	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	31kDa/46kDa
Gene Name:	sorting nexin 32
Database Link:	<a href="#">Entrez Gene 254122 Human Q86XE0</a>
Background:	May be involved in several stages of intracellular trafficking.
Synonyms:	DKFZp761P1320; FLJ30934; MGC42112; MGC57276; SNX6B




[View online »](#)

## Product images:



Western blot analysis of extracts of various cell lines, using SNX32 antibody (TA381862) at 1:1000 dilution. | Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. | Lysates/proteins: 25ug per lane. | Blocking buffer: 3% nonfat dry milk in TBST. | Detection: ECL Basic Kit. | Exposure time: 90s.



Immunofluorescence analysis of U2OS cells using SNX32 antibody (TA381862) at dilution of 1:100. Blue: DAPI for nuclear staining.