

Product datasheet for TA382119

STIP1 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: ICC/IF, IHC, IP, WB

Recommended Dilution: WB.1:500 - 1:2000

IHC,1:50 - 1:200 IF,1:50 - 1:200 IP,1:20 - 1:50

Reactivity: Human, Mouse, Rat, Monkey

Modifications: Unmodified

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of

human STIP1 (NP_006810.1).

Formulation: Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Concentration: lot specific

Purification: Affinity purification

Conjugation: Unconjugated

Storage: Store at -20°C. Avoid freeze / thaw cycles.

Stability: Shelf life: one year from despatch.

Predicted Protein Size: 59kDa/62kDa/68kDa

Gene Name: stress induced phosphoprotein 1

Database Link: Entrez Gene 10963 Human

P31948



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com




Background:

STIP1 is an adaptor protein that coordinates the functions of HSP70 (see HSPA1A; MIM 140550) and HSP90 (see HSP90AA1; MIM 140571) in protein folding. It is thought to assist in the transfer of proteins from HSP70 to HSP90 by binding both HSP90 and substrate-bound HSP70. STIP1 also stimulates the ATPase activity of HSP70 and inhibits the ATPase activity of HSP90, suggesting that it regulates both the conformations and ATPase cycles of these chaperones (Song and Masison, 2005 [PubMed 16100115]).[supplied by OMIM, Jul 2009]


Synonyms:

HOP; IEF-SSP-3521; P60; STI1; STI1L

Product images:




Western blot analysis of extracts of various cell lines, using STIP1 antibody (TA382119) at 1:1000 dilution. |Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. |Lysates/proteins: 25ug per lane. |Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence analysis of HeLa cells using STIP1 antibody (TA382119). Blue: DAPI for nuclear staining.





Immunoprecipitation analysis of 200ug extracts of HeLa cells using 1ug STIP1 antibody (TA382119). Western blot was performed from the immunoprecipitate using STIP1 antibody (TA382119) at a dilition of 1:1000.