

## **Product datasheet for TA382465**

## 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com

OriGene Technologies, Inc.

EU: info-de@origene.com CN: techsupport@origene.cn

## TGF alpha (TGFA) Rabbit Polyclonal Antibody

**Product data:** 

**Product Type:** Primary Antibodies

Applications: WE

**Recommended Dilution:** WB,1:500 - 1:2000

**Reactivity:** Human, Mouse

Modifications: Unmodified

**Host:** Rabbit

Isotype: IgG

**Clonality:** Polyclonal

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 24-98 of

human TGFA (NP\_001093161.1).

**Formulation:** Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

**Concentration:** lot specific

**Purification:** Affinity purification

Conjugation: Unconjugated

Storage: Store at -20°C. Avoid freeze / thaw cycles.

**Stability:** Shelf life: one year from despatch.

Predicted Protein Size: 16kDa/17kDa

**Gene Name:** transforming growth factor alpha

Database Link: Entrez Gene 7039 Human

P01135

**Background:** This gene encodes a growth factor that is a ligand for the epidermal growth factor receptor,

which activates a signaling pathway for cell proliferation, differentiation and development. This protein may act as either a transmembrane-bound ligand or a soluble ligand. This gene has been associated with many types of cancers, and it may also be involved in some cases of cleft lip/palate. Alternatively spliced transcript variants encoding different isoforms have been


found for this gene.

Synonyms: TFGA





## **Product images:**



Western blot analysis of extracts of various cell lines, using TGFA antibody (TA382465) at 1:500 dilution. |Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. |Lysates/proteins: 25ug per lane. |Blocking buffer: 3% nonfat dry milk in TBST. |Detection: ECL Basic Kit. |Exposure time: 180s.