

Product datasheet for **TA382595S**

TMEM106B Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB,1:1000 - 1:2000
Reactivity:	Human, Mouse, Rat
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant protein of human TMEM106B.
Formulation:	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	42kDa
Gene Name:	transmembrane protein 106B
Database Link:	Entrez Gene 54664 Human Q9NUM4
Background:	TMEM106B (transmembrane protein 106B) is a protein-coding gene. An important paralog of this gene is TMEM106A.
Synonyms:	FLJ11273; MGC33727; OTTHUMP00000155197



[View online »](#)

Product images:



Western blot analysis of extracts of various cell lines, using TMEM106B antibody ([TA382595]) at 1:1000 dilution. | Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. | Lysates/proteins: 25ug per lane. | Blocking buffer: 3% nonfat dry milk in TBST. | Detection: ECL Basic Kit . | Exposure time: 5s.