

Product datasheet for **TA382921S**

TTC5 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	WB, 1:500 - 1:2000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Reactivity:	Human, Mouse, Rat
Modifications:	Unmodified
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	49kDa
Gene Name:	tetratricopeptide repeat domain 5
Database Link:	Entrez Gene 91875 Human Q8N0Z6
Background:	Predicted to enable DNA binding activity and chromatin binding activity. Predicted to be involved in DNA repair. Predicted to act upstream of or within positive regulation of transcription by RNA polymerase II. Located in nucleoplasm.
Synonyms:	Strap



[View online »](#)

Product images:



Western blot analysis of various lysates using TTC5 Rabbit pAb (TA382921S) at 1:3000 dilution. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 90s.