

## Product datasheet for **TA382993**

### UBB Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	Dot
Recommended Dilution:	WB,1:500 - 1:2000 DB,1:500 - 1:2000
Reactivity:	Human, Mouse, Rat, Other (Wide Range)
Modifications:	K63-linkage Specific Polyubiquitin
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	A synthetic peptide corresponding to a sequence containing polyubiquitinated protein (K63 linkage).
Formulation:	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.



[View online »](#)

## Product images:



Dot-blot analysis of all sorts of peptides using K63-linkage Specific Polyubiquitin antibody (TA382993) at 1:1000 dilution.