

Product datasheet for **TA383596S**


CBL2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	WB, 1:500 - 1:2000 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Reactivity:	Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration:	lot specific
Purification:	Affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid freeze / thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	49kDa
Gene Name:	zinc finger protein 645
Database Link:	Q8N7E2
Background:	This gene encodes a member of the zinc finger domain-containing protein family. This family member contains both a RING-type and a C2H2-type of zinc finger domain, and it may function as an E3 ubiquitin-protein ligase. Protein localization suggests a role in human sperm production and quality control.
Synonyms:	FLJ25735; HAKAIL


[View online »](#)

Product images:



Western blot analysis of various lysates using ZNF645 Rabbit pAb ([TA383596]) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 180s.