

## Product datasheet for **TA383608S**


### ZNF766 Rabbit Polyclonal Antibody

#### Product data:

|                         |                                                                                                                                                                      |
|-------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:           | Primary Antibodies                                                                                                                                                   |
| Applications:           | ELISA, WB                                                                                                                                                            |
| Recommended Dilution:   | WB, 1:500 - 1:2000<br>ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.             |
| Reactivity:             | Mouse, Rat                                                                                                                                                           |
| Modifications:          | Unmodified                                                                                                                                                           |
| Host:                   | Rabbit                                                                                                                                                               |
| Isotype:                | IgG                                                                                                                                                                  |
| Clonality:              | Polyclonal                                                                                                                                                           |
| Formulation:            | Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.                                                                                                           |
| Concentration:          | lot specific                                                                                                                                                         |
| Purification:           | Affinity purification                                                                                                                                                |
| Conjugation:            | Unconjugated                                                                                                                                                         |
| Storage:                | Store at -20°C. Avoid freeze / thaw cycles.                                                                                                                          |
| Stability:              | Shelf life: one year from despatch.                                                                                                                                  |
| Predicted Protein Size: | 55kDa                                                                                                                                                                |
| Gene Name:              | zinc finger protein 766                                                                                                                                              |
| Database Link:          | <a href="#">Q5HY98</a>                                                                                                                                               |
| Background:             | Predicted to enable DNA binding activity and metal ion binding activity. Predicted to be involved in regulation of transcription, DNA-templated. Located in nucleus. |
| Synonyms:               | ZNF766                                                                                                                                                               |

[View online »](#)

## Product images:



Western blot analysis of various lysates using ZNF766 Rabbit pAb ([TA383608]) at 1000 dilution. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit . Exposure time: 1s.