

Product datasheet for **TA383704**

ALDH1L1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1/500-1/2000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	A synthesized peptide derived from human ALDH1L1
Formulation:	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	99kDa
Gene Name:	aldehyde dehydrogenase 1 family member L1
Database Link:	Entrez Gene 10840 Human O75891
Background:	Swiss-Prot Acc.O75891.Highly expressed in liver, pancreas and kidney.
Synonyms:	10-FTHFDH; DKFZp781N0997; FTHFD



[View online »](#)

Product images:



Western blot analysis of ALDH1L1 in 293T lysates using ALDH1L1 antibody.