

## Product datasheet for **TA383729S**

### Annexin V (ANXA5) Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1/5000
Reactivity:	Human, Monkey
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Purified recombinant human Annexin V protein fragments expressed in E.coli.
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% Proclin 300, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	HRP
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	30kDa
Gene Name:	annexin A5
Database Link:	<a href="#">Entrez Gene 308 Human P08758</a>
Background:	Swiss-Prot Acc.P08758.This protein is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex,which is involved in the blood coagulation cascade.
Synonyms:	Annexin-5; ANX5; CBP-I; ENX2; PAP-I; PP4; VAC-alpha



[View online »](#)

## Product images:



Western blot detection of Annexin V in HeLa and COS7 cell lysates using Annexin V mouse mAb conjugated with HRP (1:5000-1:100000 diluted), incubation time: 25u2103 1h. Predicted band size: 30KDa. Observed band size: 30KDa.