

## Product datasheet for **TA383815**

### **bcl 6 (BCL6) Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	IF: 1/50-200 WB: 1/500-1/2000 IHC-p: 1/100-1/300 ELISA: 1/20000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human BCL6. AA range:271-320
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Observed MW (kDa):78
Gene Name:	B-cell CLL/lymphoma 6
Database Link:	<a href="#">Entrez Gene 604 Human P41182</a>
Background:	Swiss-Prot Acc.P41182.
Synonyms:	BCL-5; BCL-6; BCL5; BCL6A; LAZ3; ZBTB27; ZNF51



[View online »](#)

## Product images:



Western blot analysis of Bcl6 in 293 lysates using Bcl6 antibody.