

## Product datasheet for **TA383845**

### CACNG3 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 1/1000-2000 IHC: 1/200-500
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic Peptide of CACNG3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Observed MW (kDa):30
Gene Name:	calcium voltage-gated channel auxiliary subunit gamma 3
Database Link:	<a href="#">Entrez Gene 10368 Human O60359</a>
Background:	Swiss-Prot Acc.O60359.
Synonyms:	Cacng2



[View online »](#)

## Product images:



Western blot analysis of CACNG3 in Human Brain, mouse Brain, rat Brain lysates using CACNG3 antibody.