

Product datasheet for **TA383846**

CACNG4 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 1/1000-2000 IHC: 1/200-500
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic Peptide of L-type Ca ⁺⁺ CP γ 4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Observed MW (kDa):35
Gene Name:	calcium voltage-gated channel auxiliary subunit gamma 4
Database Link:	Entrez Gene 27092 Human Q9NUBN1
Background:	Swiss-Prot Acc.Q9NUBN1.
Synonyms:	MGC11138; MGC24983



[View online »](#)

Product images:



Western blot analysis of Cav gamma 4 in Human Brain, mouse Brain, rat Brain lysates using Cav gamma 4 antibody.