

## Product datasheet for **TA383860**

### **CAMK2A Rabbit Monoclonal Antibody [Clone ID: R05-9G5]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	R05-9G5
Applications:	WB
Recommended Dilution:	WB: 1/1000
Reactivity:	Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthetic peptide of human CaMKII
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 54 kDa; Observed MW: 54 kDa
Gene Name:	calcium/calmodulin dependent protein kinase II alpha
Database Link:	<a href="#">Q13554</a>




[View online »](#)

**Background:**

Swiss-Prot Acc.Q13554. Calcium/calmodulin-dependent protein kinase that functions autonomously after Ca<sup>2+</sup>/calmodulin-binding and autophosphorylation, and is involved in dendritic spine and synapse formation, neuronal plasticity and regulation of sarcoplasmic reticulum Ca<sup>2+</sup> transport in skeletal muscle. In neurons, plays an essential structural role in the reorganization of the actin cytoskeleton during plasticity by binding and bundling actin filaments in a kinase-independent manner. This structural function is required for correct targeting of CaMK2A, which acts downstream of NMDAR to promote dendritic spine and synapse formation and maintain synaptic plasticity which enables long-term potentiation (LTP) and hippocampus-dependent learning. In developing hippocampal neurons, promotes arborization of the dendritic tree and in mature neurons, promotes dendritic remodeling. Participates in the modulation of skeletal muscle function in response to exercise. In slow-twitch muscles, is involved in regulation of sarcoplasmic reticulum (SR) Ca<sup>2+</sup> transport and in fast-twitch muscle participates in the control of Ca<sup>2+</sup> release from the SR through phosphorylation of triadin, a ryanodine receptor-coupling factor, and phospholamban (PLN/PLB), an endogenous inhibitor of SERCA2A/ATP2A2.

**Synonyms:**

CAMKA; CaMKIINalpha; KIAA0968

**Product images:**

Western blot analysis of CaMKII in mouse brain lysates using CaMKII beta antibody.