

Product datasheet for **TA383884S**

Casp9 Rabbit Monoclonal Antibody [Clone ID: R04-8I4]


Product data:

Product Type:	Primary Antibodies
Clone Name:	R04-8I4
Applications:	IP, WB
Recommended Dilution:	WB: 1/1000 IP: 1/20
Reactivity:	Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein of mouse Caspase-9
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 50 kDa; Observed MW: 50 kDa
Gene Name:	caspase 9
Database Link:	Entrez Gene 12371 Mouse Q8C3Q9
Background:	Swiss-Prot Acc.Q8C3Q9.Caspase-9
Synonyms:	APAF-3; APAF3; CASP-9; CASPASE-9c; ICE-LAP6; MCH6




[View online »](#)

Product images:



Western blot analysis of Caspase9 in 3T3 lysates using Caspase 9 antibody.



Western blot detection of Caspase-9 in C6 cell lysates using Caspase-9 Rabbit mAb(1:1000 diluted).Predicted band size:50kDa.Observed band size:50kDa.