

Product datasheet for **TA383944**


CD79A Mouse Monoclonal Antibody [Clone ID: 10A10-10C3-1H9]

Product data:

Product Type:	Primary Antibodies
Clone Name:	10A10-10C3-1H9
Applications:	IHC
Recommended Dilution:	IHC-p: 1/100-1/500
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a, kappa
Clonality:	Monoclonal
Immunogen:	Synthesized peptide derived from human CD79a
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Gene Name:	CD79a molecule
Database Link:	Entrez Gene 973 Human P11912
Background:	Swiss-Prot Acc.P11912.
Synonyms:	Ig-alpha; IGA; MB-1; MB1



[View online »](#)

Product images:

Immunohistochemical analysis of paraffin-embedded Tonsil. 1, Antibody was diluted at 1:200(4°,overnight). 2, Citrate buffer of pH6.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).