

Product datasheet for TA383976

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

CAF1 (CHAF1A) Rabbit Monoclonal Antibody [Clone ID: R06-3A5]

Product data:

Product Type: Primary Antibodies

Clone Name: R06-3A5

Applications: IF, IHC, IP, WB **Recommended Dilution:** WB: 1/1000

IHC: 1/50

ICC/IF: 1/200-1/1000

IP: 1/20

Reactivity: Human
Host: Rabbit
Isotype: IgG

Clonality: Monoclonal

Immunogen: A synthetic peptide of human p150 CAF1

Formulation: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA

Concentration: lot specific

Purification: Affinity Purified
Conjugation: Unconjugated

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Stability: 1 year

Predicted Protein Size: Calculated MW: 107 kDa; Observed MW: 150 kDa

Gene Name: chromatin assembly factor 1 subunit A

Database Link: Entrez Gene 10036 Human

Q13111






Background:

Swiss-Prot Acc.Q13111.Core component of the CAF-1 complex, a complex thought to mediate chromatin assembly in DNA replication and DNA repair. Assembles histone octamers onto replicating DNA in vitro. CAF-1 performs the first step of the nucleosome assembly process, bringing newly synthesized histones H3 and H4 to replicating DNA; histones H2A/H2B can bind to this chromatin precursor subsequent to DNA replication to complete the histone octamer. CHAF1A binds to histones H3 and H4. It may play a role in heterochromatin maintenance in proliferating cells by bringing newly synthesized cbx proteins to heterochromatic DNA replication foci .


Synonyms:

CAF; CAF-1; CAF-Ip150; CAF1; CAF1B; CAF1P150; hp150; MGC71229; P150

Product images:




Western blot analysis of p150 CAF1 in K562 lysates using p150 CAF1 antibody.



Immunohistochemistry analysis of paraffinembedded Human lung cancer using p150 CAF1/CAF antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.





Immunocytochemistry analysis of p150 CAF1(green) in Hela using p150 CAF1 antibody,and DAPI(blue)