

Product datasheet for **TA384085M**

DAZAP1 Rabbit Monoclonal Antibody [Clone ID: R05-8J7]


Product data:

Product Type:	Primary Antibodies
Clone Name:	R05-8J7
Applications:	IP, WB
Recommended Dilution:	WB: 1/2000-1/10000 IP: 1/20
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein of human DAZAP1
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 43 kDa; Observed MW: 43 kDa
Gene Name:	DAZ associated protein 1
Database Link:	Entrez Gene 26528 Human Q96EP5
Background:	Swiss-Prot Acc.Q96EP5.DAZ-associated protein 1
Synonyms:	MGC19907



[View online »](#)

Product images:



Western blot detection of DAZAP1 in K562,C6,3T3,HeLa cell lysates using DAZAP1 Rabbit mAb(1:1000 diluted).Predicted band size:43kDa.Observed band size:43kDa.