

## Product datasheet for **TA384096S**

### ICAD (DFFA) Rabbit Monoclonal Antibody [Clone ID: R09-4J6]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	R09-4J6
Applications:	WB
Recommended Dilution:	WB: 1/2000-1/10000
Reactivity:	Human, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthetic peptide of human ICAD
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 37 kDa; Observed MW: 45 kDa
Gene Name:	DNA fragmentation factor subunit alpha
Database Link:	<a href="#">Entrez Gene 1676 Human</a> <a href="#">O00273</a>
Background:	Swiss-Prot Acc.O00273.Inhibitor of the caspase-activated DNase (DFF40).
Synonyms:	DFF-45; DFF1; DFF45; ICAD



[View online »](#)

## Product images:



Western blot analysis of ICAD in rat Brain, MCF-7, C6 lysates using ICAD antibody.