

Product datasheet for **TA384099S**

DGKZ Rabbit Monoclonal Antibody [Clone ID: R04-1H8]


Product data:

Product Type:	Primary Antibodies
Clone Name:	R04-1H8
Applications:	WB
Recommended Dilution:	WB: 1/1000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein of human DGKZ/DGK-zeta
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW:104 kDa,Observed MW:124,114 kDa
Gene Name:	diacylglycerol kinase zeta
Database Link:	Entrez Gene 8525 Human Q13574
Synonyms:	DAGK5; DAGK6; DGK-zeta; hDGKzeta



[View online »](#)

Product images:



Western blot analysis of DGKZ/DGKzeta in HeLa, A549, HL-60, U2OS lysates using DGKZ antibody.