

Product datasheet for **TA384100M**


DGKZ Rabbit Monoclonal Antibody [Clone ID: R04-2A0]

Product data:


Product Type:	Primary Antibodies
Clone Name:	R04-2A0
Applications:	IP, WB
Recommended Dilution:	WB: 1/1000 IP: 1/20
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein of human DGKZ/DGK-zeta
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW:104 kDa,Observed MW:124,114 kDa
Gene Name:	diacylglycerol kinase zeta
Database Link:	Entrez Gene 8525 Human Q13574
Synonyms:	DAGK5; DAGK6; DGK-zeta; hDGKzeta


[View online »](#)

Product images:



Western blot analysis of DGKZ/DGKzeta in mouse brain, rat brain lysates using DGKZ antibody.



Western blot analysis of DGKZ/DGKzeta in Hela, A549, HL-60, U2OS, U87-MG lysates using DGKZ/DGKzeta antibody.