

Product datasheet for **TA384101**

DHRS1 Rabbit Monoclonal Antibody [Clone ID: R01-9E4]


Product data:

| | |
|-------------------------|--|
| Product Type: | Primary Antibodies |
| Clone Name: | R01-9E4 |
| Applications: | WB |
| Recommended Dilution: | WB: 1/1000 |
| Reactivity: | Rat |
| Host: | Rabbit |
| Isotype: | IgG |
| Clonality: | Monoclonal |
| Immunogen: | Recombinant protein of human DHRS1 |
| Formulation: | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA |
| Concentration: | lot specific |
| Purification: | Affinity Purified |
| Conjugation: | Unconjugated |
| Storage: | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| Stability: | 1 year |
| Predicted Protein Size: | Calculated MW: 34 kDa; Observed MW: 34 kDa |
| Gene Name: | dehydrogenase/reductase 1 |
| Database Link: | Q96LJ7 |
| Background: | Swiss-Prot Acc.Q96LJ7.SitesFeature keyPosition(s)Description |
| Synonyms: | DKFZp586l0523; FLJ14250; FLJ25430; MGC20204; SDR19C1 |



[View online »](#)

Product images:



Western blot analysis of DHRS1 in C6 lysates using DHRS1 antibody.