

Product datasheet for **TA384101M**


DHRS1 Rabbit Monoclonal Antibody [Clone ID: R01-9E4]

Product data:

Product Type:	Primary Antibodies
Clone Name:	R01-9E4
Applications:	WB
Recommended Dilution:	WB: 1/1000
Reactivity:	Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein of human DHRS1
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 34 kDa; Observed MW: 34 kDa
Gene Name:	dehydrogenase/reductase 1
Database Link:	Q96LJ7
Background:	Swiss-Prot Acc.Q96LJ7.SitesFeature keyPosition(s)Description
Synonyms:	DKFZp586l0523; FLJ14250; FLJ25430; MGC20204; SDR19C1


[View online »](#)

Product images:



Western blot analysis of DHRS1 in C6 lysates using DHRS1 antibody.