

Product datasheet for **TA384105**

digoxin Mouse Monoclonal Antibody [Clone ID: 10C5-1C5-2E7]


Product data:

Product Type:	Primary Antibodies
Clone Name:	10C5-1C5-2E7
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 1/10000
Reactivity:	Broad
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Digoxin
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% Proclin 300, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year



[View online »](#)

Product images:



Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain Tissue with EAAT2 Rabbit pAb diluted at 1:2,000.