

Product datasheet for **TA384114**

DNMT1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	FC, WB
Recommended Dilution:	WB: 1/1000 FC: 1/100
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	A synthetic peptide of human Dnmt1
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 183 kDa; Observed MW: 183 kDa
Gene Name:	DNA (cytosine-5-)-methyltransferase 1
Database Link:	Entrez Gene 1786 Human P26358
Background:	Swiss-Prot Acc.P26358.
Synonyms:	AIM; CXXC9; DNMT; FLJ16293; M.Hsal; MCMT; MGC104992



[View online »](#)

Product images:



Western blot detection of DNMT1 in Jurkat cell lysates using DNMT1 Rabbit pAb(1:1000 diluted). Predicted band size:183KDa.Observed band size:183KDa.