

Product datasheet for **TA384207M**


FANCB Rabbit Monoclonal Antibody [Clone ID: R05-8B7]

Product data:

Product Type:	Primary Antibodies
Clone Name:	R05-8B7
Applications:	WB
Recommended Dilution:	WB: 1/1000
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthetic peptide of human FANCB
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 98 kDa; Observed MW: 98 kDa
Gene Name:	Fanconi anemia complementation group B
Database Link:	Entrez Gene 2187 Human Q8NB91
Background:	Swiss-Prot Acc.Q8NB91.DNA repair protein required for FANCD2 ubiquitination.
Synonyms:	FA2; FAAP90; FAAP95; FAB; FACB; FLJ34064

[View online »](#)

Product images:



Western blot detection of FANCB in 3T3, HeLa cell lysates using FANCB Rabbit mAb (1:1000 diluted). Predicted band size: 98 kDa. Observed band size: 98 kDa.