

## Product datasheet for **TA384234**


### Flotillin 1 (FLOT1) Mouse Monoclonal Antibody [Clone ID: 6H9-1D9-3E10]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	6H9-1D9-3E10
Applications:	WB
Recommended Dilution:	WB: 1/1000-2000
Reactivity:	Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Synthetic Peptide of Flotillin-1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% Proclin 300, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Observed MW (kDa):49
Gene Name:	flotillin 1
Database Link:	<a href="#">O75955</a>
Background:	Swiss-Prot Acc.O75955.May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles.
Synonyms:	FLOT1



[View online »](#)

**Product images:**


Western blot detection of Flotillin-1 in Rat Brain, Mouse Brain cell lysates using Flotillin-1 (6H9) Mouse mAb(1:1000 diluted). Predicted band size: 49KDa. Observed band size: 49KDa.



Western blot detection of Flotillin-1 (6H9) in 3T3, Mouse brain, Rat brain, Rat hippocampus, A431 cell lysates using Flotillin-1 (6H9) antibody. Observed band size ≈ 49kDa.