

## Product datasheet for TA384268M


### FUBP1 Rabbit Monoclonal Antibody [Clone ID: R01-5D4]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	R01-5D4
Applications:	IF, IHC, IP, WB
Recommended Dilution:	WB: 1/1000 IHC: 1/500 ICC/IF: 1/50 IP: 1/20
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthetic peptide of human FUBP1
Formulation:	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	Calculated MW: 68 kDa; Observed MW: 74 kDa
Gene Name:	far upstream element binding protein 1
Database Link:	<a href="#">Entrez Gene 8880 Human Q96AE4</a>
Background:	Swiss-Prot Acc.Q96AE4.Regulates MYC expression by binding to a single-stranded far-upstream element (FUSE) upstream of the MYC promoter. May act both as activator and repressor of transcription.
Synonyms:	FBP; FUBP


[View online »](#)


## Product images:




Western blot detection of FUBP1/GBP in Mouse spleen lysates using FUBP1/GBP antibody. Predicted band size: 68kDa. Observed band size: 74kDa.



Western blot detection of FUBP1/GBP in K562, Rat Brain, HeLa cell lysates using FUBP1/GBP Rabbit mAb (1:1000 diluted). Predicted band size: 68kDa. Observed band size: 74kDa.



Immunohistochemistry of FUBP1 in paraffin-embedded Human colon cancer tissue using FUBP1 Rabbit mAb at dilution 1/100



Immunocytochemistry analysis of FUBP1 (green) in 293 using FUBP1 antibody, and DAPI (blue).