

Product datasheet for **TA384271S**

FYB1 Mouse Monoclonal Antibody [Clone ID: 2H9-H6-G10]


Product data:

Product Type:	Primary Antibodies
Clone Name:	2H9-H6-G10
Applications:	WB
Recommended Dilution:	WB: 1/500
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Purified recombinant human FYB protein fragments expressed in E.coli.
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% Proclin 300, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	120kDa
Gene Name:	FYN binding protein
Database Link:	Entrez Gene 2533 Human O15117
Background:	Swiss-Prot Acc.O15117.
Synonyms:	ADAP; FYB-120/130; p120/p130; PRO0823; SLAP-130; SLAP130



[View online »](#)

Product images:



Western blot analysis of FYB in 10, 15 and 20ug Jurkat lysates using FYB antibody.