

Product datasheet for **TA384276S**

GABARAP Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1/1000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide.
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.03% Proclin 300, pH 7.3.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Stability:	1 year
Predicted Protein Size:	14/16kDa
Gene Name:	GABA type A receptor-associated protein
Database Link:	Entrez Gene 11337 Human O95166
Background:	Swiss-Prot Acc.O95166.May play a role in intracellular transport of GABA(A) receptors and its interaction with the cytoskeleton.
Synonyms:	FLC3B; FLJ25768; MGC120154; MGC120155; MM46




[View online »](#)

Product images:



Western blot analysis of GABARAP in A549, 3T3 and HeLa lysates using GABARAP antibody.



Western blot detection of GABARAP in A549, 3T3 and HeLa cell lysates using GABARAP Rabbit pAb (1:1000 diluted). Predicted band size: 14 kDa. Observed band size: 14, 16 kDa.